




Typical Pulses for the 5mm CryoProbe™ 700MHz and 800MHz

1.	Purpose	3
2.	Scope	3
3.	Reference to Documents.....	3
4.	Implementation.....	3
4.1	Introduction.....	3
4.2	Pulse and power recommendations	4
4.3	Notes.....	5



History:

<i>Index</i>	<i>Date</i>	<i>Changes</i>	<i>Visum</i>
00	21.10.2004	document created	DAM/DMO



1. Purpose

Safe operation of a NMR probe requires that the power and the durations of RF pulses do not exceed specified values. This document provides values for safe operation of the 5mm TXI CryoProbe™ at 700MHz and 800MHz.

2. Scope

All 5mm TXI CryoProbe™ at 700MHz and 800MHz.

3. Reference to Documents

None.

4. Implementation

4.1 Introduction

The following table shows values for pulses, which provide safe operation¹⁾ with a Bruker 5mm CryoProbe™.

The given values are related to the power handling of a 5mm CryoProbe™. The effect of possible heating of the sample is not taken into account.

Please note that the given values provide safe operation only, when

- standard pulse programs from the Bruker pulse program library are used, and
- The repetition rate, which is about **d1+aq** has a minimum period of 1s.

Certain segments in a pulse sequence, like a spin lock or a decoupling sequence, are applied at a given RF power level. The maximum allowed power for those segments is expressed by the corresponding pulse length. Example:

¹H TOCSY spin lock: 200 msec @ 25 mW. This means that the maximum length of the spin lock sequence is 200 msec, and that the allowed ¹H power level corresponds to a 90 deg. pulse width of 25 mW. A power level which corresponds to a 90 deg. pulse which is shorter than 25 mW is not allowed.



4.2 Pulse and power recommendations

Note: All values require a minimum repetition period $d1+aq$ of 1sec ²⁾

Important: observe the separate notes on the following pages!

	5mm 700MHz TXI CryoProbe™	5mm 800MHz TXI CryoProbe™
¹ H		
hard pulse ³⁾ (max. length 360°)	8.0 μsec	8.0 μsec
hard pulse for lossy samples	Power level corresponding to 8.0 μsec pulse for non-lossy sample	Power level corresponding to 8.0 μsec pulse for non-lossy sample
trim pulse p28 ⁴⁾	2 msec @ 10 μsec	2 msec @ 10 μsec
TOCSY spin lock	200 msec @ 25 μsec 400 msec @ 35 μsec	200 msec @ 25 μsec 400 msec @ 35 μsec
ROESY spin lock	Up to CW for power level corresponding to a 80 μsec pulse	Up to CW for power level corresponding to a 80 μsec pulse
WALTZ16 decoupling during 13C-detection	Up to CW for power level corresponding to a 100 μsec pulse	Up to CW for power level corresponding to a 100 μsec pulse
DIPS12-decoupling in triple resonance	400 msec @ 35 μsec	400 msec @ 35 μsec
¹³ C		
hard pulse ⁶⁾ (max. length 360°)	15.0 μsec	15.0 μsec
trim pulse ⁵⁾	2msec @ 22 μsec	2msec @ 22 μsec
CC spin lock ⁵⁾	20msec @ 22 μsec	20msec @ 22 μsec
GARP decoupling ⁷⁾	140msec @ 55 μsec	140msec @ 55 μsec
selective pulses ⁸⁾	G4: 352 μsec Q3: 220 μsec CHIRP: 2ms @ 25 μsec	G4: 308 μsec Q3: 210 μsec CHIRP: 2ms @ 25 μsec
¹⁵ N		
hard pulse ⁶⁾ (max. length 360°)	40.0 μsec	40.0 μsec
GARP4 decoupling ⁷⁾	140 msec @ 170 μsec	140 msec @ 170 μsec
CPMG T2 ⁹⁾	250 msec @ 80 μsec (see warning ⁹⁾)	250 msec @ 80 μsec (see warning ⁹⁾)
² H		
hard pulse (max. length 360°)	150 μsec	150 μsec
WALTZ16 decoupling	100 msec @ 250 μsec	100 msec @ 250 μsec
Z-Gradient		
Absolute max. current ¹⁰⁾	10A	10A
Max. overall length ¹¹⁾	10ms @ 10A	10ms @ 10A



Table 1

- 1) Although the values given do not necessarily define the actual limits, these conditions have been determined to provide safe operation and are valid unless superseded in whole or in part by (1) a more recent version of this document, (2) a factory document pertaining to the specific probe, or (3) written approval by BRUKER BioSpin AG. Operation in excess of the conditions given is not permitted and will void the warranty as it may result in a damage to the CryoProbe.
- 2) Repetition period related issues:
- In any case, the power level for a given pulse length (which must be equal or shorter than the specified pulse length value) must never exceed the power maximum given in Table 1.
 - The pulse length values given in Table 1 are based on a repetition period of approximately 1sec. This defines and limits the average power to which the probe may be exposed. Operation with longer repetition periods is allowed without additional restrictions, provided the pulse lengths are still limited according to the values given in Table 1.
- 3) Non-lossy samples are assumed, such as organic solvents or water without salt.
- 4) A number of inverse correlation experiments use the proton trim pulse *p28*, like the pulse program *invietgpsi* (*XWINNMR up to version 3.1*) or *hsqcetgpsi* (*XWINNMR version 3.5*). According to our experience this trim pulse is only rarely required. Verify if any improvement is necessary or obtained using this pulse.
- 5) An experiment might contain a spin lock followed by a trim pulse. Note that the total length of spin lock and trim pulse must not exceed the value which is allowed for the spin lock only.
- 6) When hard pulses are applied simultaneously on the ^{13}C and ^{15}N channel, the power on both, the ^{13}C and the ^{15}N channel, has to be reduced by 3dB each.
- 7) Simultaneous decoupling of ^{13}C and ^{15}N requires that reduced power is used on both channels. 3dB less power should be used for both ^{13}C and ^{15}N . For simultaneous decoupling it is recommended to use optimised conditions for adiabatic decoupling on ^{13}C .
- 8) Pulses Q5 and Q3 are used in triple resonance experiments. The 2ms CHIRP pulse is used in HSQC experiments for refocusing. For shaped selective pulses there are usually no restrictions, as long as the peak power of the shape does not exceed the power limit of a hard pulse. When using non-standard pulses, in addition to the peak power, it should be also verified that the averaged power over the pulse length does not exceed any of the power limits (compare to the listed limits for pulses with a comparable duration).



⁹⁾ The pulse sequence *invit2etf3gpsi* (*XWINNMR* up to version 3.1) or *hsqct2etf3gpsi* (*XWINNMR* version 3.5). The length of the CPMG element is *d20* and should not exceed 250 msec. The 180 deg. pulse is *p30* at a power level *pl23*, the pulse *p30* should not be shorter than 160 μ sec (180 deg. pulse, applied at a power level corresponding to a 90 deg. pulse of 80 μ sec). The given value is true only, when the pulse *p30* is part of a CPMG sequence with an inter-pulse delay of 900 μ sec (expressed by *d21*2* in the pulse sequence *invit2etf3gpsi*). **Warning: Never apply the above power level and full uninterrupted duration in a spin lock sequence.**

¹⁰⁾ The absolute max. gradient current I_{max} is listed in Table 1.

¹¹⁾ It is permissible to apply I_{max} for the listed maximum overall length T_{max} . The overall length of the gradients is the sum of all gradients applied during a repetition period.

A gradient sequence which has, within a repetition period, an overall length T_G which is longer than T_{max} , is still permissible if the maximum gradient current is reduced to I_G , according to the following equation:

$$I_G = I_{max} * \sqrt{T_{max}/T_G}$$

In a more general formulation (using the values of $I_{max} = 10A$ and $T_{max} = 10msec$ for clarity), for N pulses within any repetition (one second) period, the following condition must hold (I_{Gi} , T_{Gi} denote the current and duration of the i-th gradient pulse):

$$\sum_{i=1}^N T_{Gi} I_{Gi}^2 \leq T_{max} * (I_{max})^2 = 10ms * (10A)^2 = 1000 \cdot msA^2$$

Examples:

1. For an overall gradient length of $T_G = 20ms$ within any one second period, the maximum gradient current I_G is 7.07A.
2. A pair of gradient pulses, consisting of a first pulse with a current of 10A and a length of 5msec, and a second pulse with a current of 5A and a length of 20msec is just permissible.
3. A DC current of 1A is just permissible.